# Instruction Manual

## Automatic Blood Pressure Monitor with PC LINK LEDs indicators for Arm V2.1

English

**BPC-81U**

# Contents

<b>Introduction.....</b>	<b>4</b>
<b>Warning Signs &amp; Safety Statements .....</b>	<b>4</b>
<b>Health and Blood Pressure Values .....</b>	<b>4</b>
<b>Warning .....</b>	<b>5</b>
<b>Suggestions before Measuring Blood Pressure .....</b>	<b>5</b>
<b>Device Descriptions.....</b>	<b>7</b>
<b>Description of Control Buttons .....</b>	<b>8</b>
<b>Description of LCD Display Symbols .....</b>	<b>8</b>
<b>Installing/Replacing Batteries .....</b>	<b>9</b>
<b>Program Date and Time .....</b>	<b>11</b>
<b>How to Apply the Arm Cuff .....</b>	<b>12</b>
<b>Correct Measurement Position .....</b>	<b>14</b>
<b>Measurement while sitting down .....</b>	<b>14</b>
<b>How to Take A Measurement .....</b>	<b>14</b>
<b>How to SEND MEASUREMENTS TO PC? .....</b>	<b>16</b>
<b>What is Irregular Heartbeat? .....</b>	<b>17</b>
<b>What is Arrhythmia? .....</b>	<b>17</b>
<b>How to Use the Memory Features .....</b>	<b>18</b>

<b>How to delete Single or All Memories.....</b>	<b>19</b>
<b>HealthPlus Software Installation and User Manual .....</b>	<b>20</b>
<b>Maintenance .....</b>	<b>26</b>
<b>Error Message &amp; Troubleshooting .....</b>	<b>26</b>
<b>TROUBLESHOOTING.....</b>	<b>28</b>
<b>Specifications.....</b>	<b>30</b>
<b>Warranty Information .....</b>	<b>31</b>

## Introduction

Thank you for purchasing the Medron BPC-81U Automatic Blood Pressure Monitor with LEDs for Arm.

For your safety, please be sure to read this Instruction Manual thoroughly before using the Blood Pressure Monitor. For other specific information about your own blood pressure, CONTACT YOUR PHYSICIAN.

## Warning Signs & Safety Statements

ATTENTION!!  
Read user's manual before use.


Type BF applied part.

## Health and Blood Pressure Values

Classification of Blood Pressure by the WHO:

The World Health Organization (WHO) and the Inter-national Society of Hypertension (ISH) developed the Blood Pressure Classification shown as below.

**According to the blood pressure classification by the WHO/ISH (revised in 1999)**

## Warning

- If the event of the cuff pumping up and not stopping, please open the cuff at once.
- Do not measure repeatedly for a long period of time.
- This device is not supposed to be used if your wrist has any wound or injury.
- If locks connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intra-vascular fluid systems, which case may allow air to be pumped into a blood vessel.
- Take batteries out of the housing if the unit is not to be used for a long time.
- Do not use different type, date or brand-name batteries, and this may lead to an incorrect reading.
- Leaving the batteries in the blood pressure monitor for a long period and not using the device may damage the main unit.


## Suggestions before Measuring Blood Pressure

1. It is recommend that persons has one of the diseases as heart troubles, kidney trouble, hyper-tension, diabetes,

arteriosclerosis or frustration problem of blood vessel, please consult your physician before and during using the blood pressure monitor to determine if your wrist blood pressure accurately reflects your actual blood pressure.

2. Please sit or lie down for at least 5-10 minutes before taking a measurement.
  3. Avoid eating, smoking, exercising, bathing and drinking (alcohol/caffeine/or tea) for at least 20-30 minutes before taking a measurement, and avoid talking during a measurement.
  4. Be sure the wrist cuff is at the same level as your heart.
  5. Stress raises blood pressure. Keep calms and relax yourself during a measurement.
  6. Wait 2-3 minutes before taking another blood pressure measurements.
  7. Keep the records for your physician as reference.
- Reminded that one measurement does not provide you with real blood pressure value, only long-term records will reveal the long-term trend of your blood pressure. Try to measure your blood pressure at the same time each day for consistency.


## Device Descriptions

## Description of Control Buttons

- **POWER/START**: To turn on or turn off the unit. Start inflating to measure blood pressure and pulse.
- **⌚ TIME**: To set the Date/Time.
- **MEMO/▲**: To review the measured record(s) manually / To select the digit of Date/Time.
- **SCAN/▼**: Automatically review all previously memorized record(s) / To select the digit of Date/Time / To delete memorized record(s).
- **PC SYNC**: To send data in memory to PC. Press this button twice to send data to PC.

## Description of LCD Display Symbols

## Installing/Replacing Batteries

You can use the unit by using Batteries.


Be sure the unit turned off before installing/replacing batteries, or the memories & Date/Time setting will be erased. If the batteries are removed or unplug the adapter for more than 40 seconds, the memories and the Date/Time setting will be erased. Then, you must reset the Date/Time.

Slide the battery compartment to remove it.

Place the batteries (4 “AA” Alkaline batteries) in the compartment as indicated by the polarity symbols (+ and -) marked inside. All the display symbols appear for 3 seconds and then return to standby mode.

Replace the battery compartment cover.

Note:

When  (Low Battery indicator) appears on the display, the batteries are weak. Change the batteries as soon as possible, or it may not be able to take any measurement.

And, the memories and the Date/Time setting of the unit may be lost.

We recommend you to use the ALKALINE batteries for the best performance and the longest life.

Remove the batteries if not to use the blood pressure

monitor for an extended period (i.e., 3 months). The batteries can leak chemicals that can destroy electronic parts.

Immediately replace run-down batteries with new ones.

## **Program Date and Time**

Continuously press TIME button for 2 seconds until “2006” shown on the screen and starts to flash.

Press MEMO/▲ button or SCAN/▼ button to set the number representing the Year.

(You can press and hold MEMO/▲ or SCAN/▼ buttons for about 2 seconds, the number of the Year will be displayed in the faster speed.)


Press TIME button to confirm the Year setting. Then the number representing the Date will flash.

In the same manner, set the Month, Date, Hour and Minute.

Note: After replacing batteries, always check the Date/Time is properly set or not first, then use the blood pressure monitor. If the batteries are removed for more than 40

seconds, you must reset the Date/Time.

## How to Apply the Arm Cuff



Insert the air plug of the arm cuff into the air jack of the blood pressure monitor.

The cuff should be assembled correctly when it is removed from the box. If it is not, pass the end of the cuff furthest from the tubing through the metal D-ring to form a loop. The smooth cloth should be on the inside of the cuff loop.

If the cuff is assembled correctly, the sewn hook material will be on the outside of the cuff loop and the metal D-ring will not touch your skin.

Put your left arm through the cuff loop. The bottom of the

cuff should be approximately 1/2" above the elbow. The cuff tubing should be positioned along the middle of the inside of your arm as shown.


1. Pull the cuff so that the top and bottom edges are tightened evenly around your arm.
2. When the cuff is positioned correctly, press the sewn hook material firmly against the pile side of the cuff.
3. The cuff should be wrapped so that it fits snugly and stays in place. You should be able to fit one finger between the cuff and your arm.

## **Correct Measurement Position**

It is extremely important for the measure values, which that the cuff is at the same level of the heart. Higher or deeper position will give deviations from the true value.

### **Measurement while sitting down**

1. Sit in a chair with your feet flat on the floor and rest your arm on a table so that the cuff is at the same level as your heart.
2. Relax your arm and turn your palm upward.
3. Make sure there are no kinks in the air tubing.
4. After taking deep breaths for several times before measurement, take a natural breath, and then relax your body.

### **How to Take A Measurement**

Please read and follow directions in the section

“Suggestions before Measuring Blood Pressure” and “How to Apply the Arm Cuff” before taking a measurement.

### **Selecting User 1 or 2**

Press  **POWER** button to turn on the unit. All the symbols displayed for a couple seconds then the number **0** and person icon will flash. Then you can press **▲** or **▼** to choose user 1 or 2. The unit is ready to take “measurements, when the icon stopping flash. If no specific user chosen, the default user is “1”.

Press **START** button, the cuff automatically inflates to approximately 170mmHg. (50mmHg higher than the average systolic blood pressure reading.)

1. When the cuff reaches your correct inflation level, it begins to deflate automatically with decreasing numbers shown on the display. The Heart Symbol  flashes as the onset of oscillation signals.

\*

. If during deflation the blood pressure monitor senses that higher pressure is required to measure the user's blood pressure, the device will stop deflation and pump again to

reach the workable pressure (up to the max. 280 mmHg). If the inflation continues and doesn't stop, please untie the cuff immediately or turn off the unit by pressing the POWER. The monitor exhausts the residual air when the measurement is completed. The Heart Symbol  stops flashing and your blood pressure and pulse will appear on the display. In the mean time, the unit announces the measured value of Systolic, Diastolic and Pulse. After taking the measurement, the latest measured readings will always be displayed before the power is turned off.

Note:

Wait for at least 3 minutes if you want to take another measurement. Waiting allows your blood vessels to bounce back to normal status. More rest time may be required between measurements depending on individual physiological conditions.

## **How to SEND MEASUREMENTS TO PC?**


Please read SOFTWARE INSTALL AND USER'S GUIDE section in this manual to install and run Medron HealthPlus

Software on PC.

## What is Irregular Heartbeat?

An irregular heartbeat is defined as a heartbeat that varies by 25% from the average of all heartbeat intervals while the unit is measuring the systolic and diastolic blood pressure.

If such an irregular rhythm is detected more than twice during measurement, the irregular heartbeat symbol appears on the display with measurement results.

## What is Arrhythmia?

A heartbeat is stimulated by electrical signals that cause the heart to contract.

Arrhythmia is a condition where the heartbeat rhythm is abnormal due to flaws in the bio-electrical system that drives the heartbeat. Typical symptoms are skipped heartbeats, premature contraction, an abnormally rapid (tachycardia) or slow (bradycardia) pulse. This can be caused by heart disease, aging, physical predisposition, stress, lack of sleep, fatigue etc.

Arrhythmia can only be diagnosed by a doctor through a special examination.

Whether the appearance of the irregular heartbeat symbol in the results indicates arrhythmia or not can only be determined by an examination and diagnosis by your doctor.

**Note:**

If the irregular heartbeat symbol is shown frequently, please make your doctor aware of it. Conducting self-diagnosis and treatment based on measurement results are dangerous. Be sure to follow the instructions of your doctor.

## **How to Use the Memory Features**

The blood pressure and pulse measured and displayed are automatically stored in the memory. The blood pressure monitor can store up to 80 measurements. If more than 80 measurements have been taken, the earliest measurement will be erased and replaced by the most recent measurement.

1. Review the Average Value and previously measurements:  
Press MEMO/▲ button, the unit first display the average values of Systolic, Diastolic and Pulse. (Averaged from all memorized readings.) Press MEMO/▲ button again to review the measurement from the latest to the earliest.
2. The latest measurement (M 01) appears when you press MEMO/▲ button. M 01 is the latest reading measured. The Systolic, Diastolic, Pulse and Date/Time of the measurement will be displayed. Press MEMO/▲ button to repeatedly review other memorized readings. The

readings will be displayed as M 02, M 03, M 04 ... etc. You can press SCAN/▼ button to switch to the automatic display mode.

3. Scan all previous measurements:

Press SCAN/▼ button, the unit will display all measurements one by one in one in one second interval. You can press MEMO/▲ button to switch to the single display mode.

4. The 3 colors LED warning lights:

It will indicate the status of each blood pressure measurement is Normal (green), Higher (yellow) or Abnormal High (reddish).

## How to delete Single or All Memories

1. Delete single measurement:

You can delete the single memory just after taking the measurement. Press SCAN/▼ button then “dEL” will appear and flash. Press SCAN/▼ button again to confirm the delete function.

2. Delete all measurements:

Turned off the unit first then keep pressing SCAN/▼ until “dEL ALL” displayed. Press SCAN/▼ button again and all memories will be deleted when “00 ALL” shown on the screen.

# HealthPlus Software Installation and User Manual

Please check

<http://www.medronhealth.com/products/BPC-81U/docs/latest> for latest updates on the software manual.

## 1. Installation

### 1.1 Software Installation

- Download latest version of Medron HealthPlus software from <http://www.medronhealth.com/products/healthplus/software/latest>
- Unzip it to known location on your computer
- Double click setup.exe file in unzipped folder. Please make sure computer is connected to Internet.
- It shall open Application Install dialog box
- Click Install button (as shown in the picture below)
- Follow instructions that will come up on your screen.
- Installation concludes with starting of HealthPlus application.  
HealthPlus application is also accessible by following Start-> Programs -> Medron -> HealthPlus

HealthPlus installation performs following tasks transparent to user,

- Installs.NET 4.x if computer does not have it already – This task


may take up to 12 minutes in case .NET is not installed on the computer. Please ensure that computer is connected to the internet.

- Installs driver for the device for Windows8, Windows7 and WindowsXP
- Installs application

## 2. User Operations

### 2.1 Settings

In this view, you can assign names and nicknames for two users. Simply double click on any fields and type in the name, DOB, gender info. A unique user id is automatically assigned for each user.

### 2.2 Sync

Click on Sync Button and follow onscreen directions. **Ensure that you click on Ignore duplicate entries box to ensure that entries from Medron Blood pressure unit are not duplicated, whenever**

you sync again.

## 2.3 Readings

Tabular representation of Systolic, Diastolic readings with time stamp date information.

### 2.3.1 Add a new reading

This feature allows you to add a reading manually. You can add a new reading by double clicking the last row that is empty (No data in that row). This will pre-populate all fields of this row with fixed set of values. Double click on each field to modify it.

### 2.3.2 Edit an existing reading

This feature allows you to edit an existing reading. Double click the left mouse button on the desired cell of any row (reading) to edit the data in that cell.

### 2.3.3 Deleting a reading

This feature allows you to delete any single row by pressing Delete button on the right end of the row. It also allows deletion of multiple rows by selecting multiple rows, and then pressing delete key on the keyboard.

The screenshot displays the MedRon HealthPlus application window. The title bar reads "HealthPlus". The interface includes a sidebar with navigation options: Home, Readings, Charts, My Diary, Sync, Report, Settings, and Help. The main content area is titled "Blood Pressure Readings" and contains a table with the following data:

Date	Time	Systolic	Diastolic	Pulse	Stress Level	Reading Notes	Actions
Tuesday, October 23, 2012	2:07 PM	120	80	70	Low	Taken with Medron BP	Details Delete
Tuesday, October 23, 2012	2:02 PM	120	80	70	Low	Taken with Medron BP	Details Delete

## 2.4 Charts

Graphical representation of Systolic, Diastolic and Pulse readings plotted over time. You can zoom in and zoom out using your mouse scroll bar. Scroll down to view graphical plot of pulse information

## 2.5 My Diary

My Diary allows user to document daily notes about general health, fitness or almost any thing that would be good for Physician to know during next visit.


The screenshot shows the 'HealthPlus' application window with the 'My Diary' section selected. The interface includes a table for logging diary entries with columns for Date, Note Title, and a text area for the note content.

Date	Note Title	Note
1/3/2012 12:00:00 AM	Note1	
1/23/2012 12:00:00 AM	Note2	


## 2.6 Report

This view allows user to print the report. Report includes all readings, systolic/diastolic, pulse reading charts, user name and relevant details needed for your physician. It supports “Print Report” as well as “Save to Excel (.xlsx)” features.

## 2.7 Help

Help view has all manuals including Quick start guide, software installation, instructional manual, user manual, support tickets form, FAQ.

## **Maintenance**

### **Storage**

1. Always put the device in the storage box attached when not in use or during transporting.
2. If you won't use the device for an extended period like over 3 months, remove the batteries.
3. Do not store the device under direct sunlight, high temperature, humid, dusty, or corrosive environment.
4. Do not stretch, bend or tug the wrist cuff too hard.
5. Do not wet the unit.

### **Cleaning**


1. Do not use an alcoholic-base, benzene, gasoline or solvent agent to clean the device.
2. Wipe with a moistened cloth with clean water or a mild cleaning agent and then dry the device with a soft and dry cloth.
3. The dirt on the cuff can be wiped by the moistened cloth and soap.
4. Do not flush the device and cuff.
5. Do not bend the cuff or sleeve and wrap the sleeve to the inside-out surface.
6. Do not disassemble the device or dismantle the cuff or try to repair by yourself. If any problem happens, bring the device to the distributor for further examination.
7. Do not operate the device under the severe environment of extreme temperature or humidity, or direct sunshine.

## **Error Message & Troubleshooting**

The unit is manufactured through the complete quality assurance system.

If the measurement is not made correctly, the following error will be displayed. Measure your blood pressure again correctly.

## ERROR MESSAGES

Error Indicators	Cause	How to correct
	<ul style="list-style-type: none"> <li>● Incorrect posture or talks during measurement.</li> <li>● There are any artificial interference and noise during measurement.</li> <li>● User is tense and does not relax or has severe heart beat problem.</li> </ul>	<ul style="list-style-type: none"> <li>● Cuff is not fastened well.</li> <li>● Wait at least 3 minutes before taking another measurement.</li> <li>● User has special characteristic.</li> </ul>
	<ul style="list-style-type: none"> <li>● No pulse or blood pressure sensed.</li> <li>● Incorrect posture during measurement.</li> <li>● Cuff is not wrapped properly.</li> </ul>	<ul style="list-style-type: none"> <li>● Cuff is not fastened correctly or not tight enough.</li> </ul>
	<ul style="list-style-type: none"> <li>● Insufficient or failure inflation.</li> <li>● Cuff broken.</li> </ul>	<ul style="list-style-type: none"> <li>● Cuff is not fastened well; please re-fasten the cuff.</li> <li>● Pump, or venting valve or deflation valve failure, please contact the distributor.</li> </ul>

If a trouble occurs while you are using the blood pressure monitor, check the attached Troubleshooting first. If you can not measure your blood pressure correctly even after taking the following countermeasures, the monitor may be broken. Please send the unit to the distributor for repair. Do not dismantle or try to repair the monitor by yourself.

## TROUBLESHOOTING

Condition	Cause	Correction
After the batteries are installed, no indication or symbol displayed on the LCD. The wrong indication when power is turned on.	<ul style="list-style-type: none"> <li>●Batteries are not installed properly.</li> <li>●Batteries run out the electricity.</li> </ul>	<ul style="list-style-type: none"> <li>●Check and correct the battery polarities.</li> <li>●Remove the old batteries and wait for one minute, then install the new batteries.</li> </ul>
The inflation cannot be done or the air pressure cannot rise.	Cuff is not wrapped properly.	<ul style="list-style-type: none"> <li>●Check the arm position and fasten the cuff correctly and re-measure the pressure again.</li> </ul>
	Battery low. Battery voltage is excessively low.	<ul style="list-style-type: none"> <li>●In order to keep the previous measurements &amp; date/time, ensure the unit is turned off before replacing batteries.</li> <li>●Replaced with new “AA” Alkaline batteries.</li> </ul>

Condition	Cause	Correction
The blood pressure cannot be taken, and the LCD shows error message or wrong result.	-----	<ul style="list-style-type: none"> <li>● Re-fasten the cuff.</li> <li>● Relax yourself.</li> <li>● Keep the cuff and heart at the same level during the measurement.</li> <li>● Keep silent and still in measurement.</li> <li>● If the patient has severe heart beat problem, and the blood pressure reading may not be correct.</li> </ul>
Under normal operating circumstance, the reading at home is different from that of the clinics.	-----	<ul style="list-style-type: none"> <li>● Usually people feel tenser in the hospital or clinic (white robe effect) than in their home.</li> <li>● Record the results and consult with your physician.</li> </ul>
Each measurement has various reading.		<ul style="list-style-type: none"> <li>● The variation might be due to the different environments.</li> <li>● The blood pressure is changing according to the physiological or psychological status of the human body.</li> </ul>

## Specifications

Model No.	BPC-81U
Display	LCD Digital Display
Measurement Method	Oscillometric method
Measurement Range	Pressure: 20~280mmHg Pulse: 40~190 / minute
Accuracy/ Calibration	Pressure: $\pm 3$ mmHg or 2% of reading Pulse: $\pm 5\%$ of reading
Memory capacity	80 sets memories (80x1 user)
Inflation	Automatic inflation by air pump
Deflation	Automatic liner pressure deflation valve
Pressure Release	Automatic exhaust valve
Sensor	Semiconductor pressure sensor
Power Supply	6V DC (4 "AA" alkaline batteries)
Measuring Interval	1 mmHg
Operating Temperature / Humidity	+10°C~+40°C (50 °F~104 °F) 40%~85% R.H. maximum
Storage Temperature / Humidity	-20°C~+50°C (-4 °F~122 °F) 15%~90% R.H. maximum
Battery Life (Alkaline batteries)	Approx. 250 uses (It is estimated on the premise that backlight hasn't been enabled.)
Unit Weight	Approx. 400grams (not including batteries)
Cuff Circumference	Fits arm circumference 8.7" ~ 12.6"
Unit dimensions	Approx. 115 (l) x 125 (w) x 50 (h) mm
Accessories	Storage Case, Instruction Manual, 4 "AA" alkaline batteries

Specifications are subjected to change.

## Warranty Information

Your product MEDRON BPC-81U, Blood Pressure Monitor, is guaranteed to be free from defects in workmanship and materials under normal use for a period of 4 years from the date of purchase.

For repair within four year warranty, our Authorized Service Agent must be advised the defects of the device. This warranty covers parts and labor only under normal usage, any transportation cost and freight incurred will be at the owner's expenses. Any defect from natural causes, e.g., flood, hurricane etc., is not within this warranty. This warranty does not cover damage incurred by: usage NOT following the instructions, accidental damage, or being tampered with or serviced by unauthorized Service Agents.

Monitor subjected to misuse, abuse, negligence of the instruction manual, unauthorized repair or modifications will waive the warranty. The warranty excludes expendables and consumables like batteries. All warranty claims must be directed to the distributor for the sale of the device. The content of this warranty is subject to change without further notice.

## BLOOD PRESSURE CHART

### High Blood Pressure Range

Systolic pressure (mm Hg)	Diastolic pressure (mm Hg)	Stages of High Blood Pressure
210	120	Stage 4
180	110	Stage 3
160	100	Stage 2
140	90	Stage 1

### Normal Blood Pressure Range

Systolic pressure (mm Hg)	Diastolic pressure (mm Hg)	Pressure Range
130	85	High Normal Blood Pressure
120	80	Normal Blood Pressure
110	75	Low Normal Blood Pressure

### Low Blood Pressure Range

Systolic pressure (mm Hg)	Diastolic pressure (mm Hg)	Pressure Range
90	60	Borderline Low blood Pressure
60	40	Too Low Blood Pressure
50	33	Dangerously Low Blood Pressure

Source for Tables on Page 32 & 33: Disabled-world.com

## BLOOD PRESSURE RANGE BY AGE

<b>Age 15 to 19</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	105	117	120
Diastolic Range	73	77	81
<b>Age 20 to 24</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	108	120	132
Diastolic Range	75	79	83
<b>Age 25 to 29</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	109	121	133
Diastolic Range	76	80	84
<b>Age 30 to 34</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	110	122	134
Diastolic Range	77	81	85
<b>Age 35 to 39</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	111	123	135
Diastolic Range	78	82	86
<b>Age 40 to 44</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	112	125	137
Diastolic Range	79	83	87
<b>Age 45 to 49</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	115	127	139
Diastolic Range	80	84	88
<b>Age 50 to 54</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	116	129	142
Diastolic Range	81	85	89
<b>Age 55 to 59</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	118	131	144
Diastolic Range	82	86	90
<b>Age 60 to 64</b>	<b>Min</b>	<b>Average</b>	<b>Max</b>
Systolic Range	121	134	147
Diastolic Range	83	87	91

## PERSONAL BLOOD PRESSURE'S RECORD

Name \_\_\_\_\_ Age \_\_\_\_\_ Weight \_\_\_\_\_

Date	Time	SYS./DIA./PULSE	Remark

## Customer Complain Method

If you have any problems with the device, please notify the distributor within the warranty period, with the:

- (1) Model no.
- (2) Serial no.
- (3) The actual problem
- (4) Quantity of the devices

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Model No.	BPC-81U
Display	LCD Digital Display
Measurement Method	Oscillometric method
Measurement Range	Pressure: 20~280mmHg Pulse: 40~190 / minute
Accuracy/ Calibration	Pressure: $\pm 3$ mmHg or 2% of reading Pulse: $\pm 5\%$ of reading
Memory capacity	80 sets memories (80x1 user)
Inflation	Automatic inflation by air pump
Deflation	Automatic liner pressure deflation valve
Pressure Release	Automatic exhaust valve
Sensor	Semiconductor pressure sensor
Power Supply	6V DC (4 "AA" alkaline batteries)
Measuring Interval	1 mmHg
Operating Temperature /Humidity	+10°C~+40°C (50 °F~104 °F) 40%~85% R.H. maximum
Storage Temperature / Humidity	-20°C~+50°C (-4 °F~122 °F) 15%~90% R.H. maximum
Battery Life (Alkaline batteries)	Approx. 250 uses (It is estimated on the premise that backlight is not enabled.)
Unit Weight	Approx. 400grams (not including bteries)
Cuff Circumference	Fits arm circumference 8.7" ~ 12.6"
Unit dimensions	Approx. 115 (l) x 125 (w) x 50 (h) mm
Accessories	Storage Case, Instruction Manual, 4 "AA" alkaline batteries

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